

## *Isolation of RNA*

- 1.) Weigh sample to between 50 and 100 mg. Don't let sample defrost!
- 2.) Put sample in motor and pestal
- 3.) Add liquid nitrogen
- 4.) Grind sample
- 5.) Pour sample and liquid nitrogen into eppendorf tube
- 6.) Wait for liquid nitrogen to evaporate
- 7.) Add **450 uL of RTL buffer** and **4.5 uL B-ME**  
(if making buffer solution ahead of time: 1 ml buffer x 10 uL B-Me)
- 8.) Vortex
- 9.) Incubate at 56 Celsius for 1-3 minutes
- 10.) Pipet all of the lysate (approximately 550 ul) into the **lilac QIA shredder** and put that into a 2 ml eppendorf tube. (For this step cut off end of pipet tip with ethanol-steralized scissors)
- 11.) **Centrifuge for 2 min** max speed. (may need to snip the pipet tip)
- 12.) Pipet the flow through supernatant to a new microcentrifuge tube without disturbing the pellet.
- 13.) Add **225 ul ethanol** to the supernatant, and mix by pipetting. Precipitate may form. No centrifuging!!! Go immediately to next step!
- 14.) Add **650 ul** of above to an **Rneasy mini column** (pink) placed in a 2 ml collection tube (supplied).
- 15.) **Centrifuge for 15 s** at > 10,000 rpm/min.
- 16.) **Discard** whatever lands in the centrifuge tube, but do not throw away the tube.
- 17.) Add **700 ul buffer RW1** to the pink column.
- 18.) **Centrifuge for 15s** at 10000 rpm. Now throw away the tube and the flow through.
- 19.) **Transfer** the pink column to a new 2 ml collection tube.
- 20.) Pipet **500 ul of buffer RPE** to the pink column.
- 21.) **Centrifuge for 15s** at 10000 rpm.
- 22.) **Discard** the flow through.
- 23.) Add **500 ul RPE** to the pink column.
- 24.) **Centrifuge 2 min** at 10000 rpm.
- 25.) **Transfer** RNeasy column into a new 2 ml tube and **centrifuge for 1 min.**
- 26.) **Repeat.**
- 27.) **Transfer** the RNeasy column to a new 1.5 ml collection tube (supplied)  
and
- 28.) Pipette **30-50 ul RNase free water** directly into the RNeasy membrane. Don't let it touch the sides.
- 28.) **Centrifuge for 1 min** at 10000 rpm.
- 29.) **Store** the RNA in a -20 - -70 celsius freezer.

### *Quantitation of RNA*

- 1.) Wash cuvette with .1M NaOH, 1mM EDTA, followed by washing with RNase-free water. ( Follow the solutions guide). It has to be stored for 12 hrs.
- 2.) 10 ul of RNA sample + 490 ul distilled water (1/50 dilution).
- 3.) Measure absorbance in a 1 ml cuvette at 260 nm.

Concentration of RNA sample= 40 x Absorbance x 50.

Run a 1 % agarose gel (If you want to). 10 ul sample + 2 ul 10x loading dye.